

## Measurement of pRSET-A DNA Concentration

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Wednesday, May 31, 2023

**Protocol:**

1. Cleaned NanoDrop spectrophotometer with DI water
2. Blanked NanoDrop with 2  $\mu$ L Elution Buffer
3. Cleaned NanoDrop before loading 2  $\mu$ L pREST vector

**Results:** The DNA concentration of the pREST vector was measured to be 5.6 ng/ $\mu$ L.

**Conclusion:** We will transform E. coli with pRSET-A in order

mScarlet is used for defluorinating enzymes

pRSET – amp resistance can't have kan; for beta oxidation genes

## **Transformation of E. coli with pRSET**

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Wednesday, May 31, 2023

**Protocol:**

1. Thawed NEB 10-beta Competent E. coli cells on ice for 10 minutes
2. Added 25  $\mu$ L E. coli cells to a microfuge tube
3. Added 2  $\mu$ L pRSET to the tube
4. Flicked tube to mix and briefly centrifuged the samples down
5. Placed tube on ice for 30 minutes
6. Heat shocked E. coli in a 42 °C water bath for 30 seconds
7. Placed tubes on ice for 5 minutes
8. Added 475  $\mu$ L outgrowth medium to each tube
9. Incubated E. coli at 37 °C and shook at 400 rpm for 1 hour

**Results:** N/A

**Conclusion:** N/A

**Did transformation twice because for the first time, they used LB-amp instead of outgrowth medium; redid**

## **Transformation of E. coli with pREST-iGluSnFR**

**You used this because pRSET A wasn't a high enough transformation ???**

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Wednesday, May 31, 2023

**Protocol:**

1. Thawed NEB 10-beta Competent E. coli cells on ice for 10 minutes
2. Added 25  $\mu$ L E. coli cells to a microfuge tube
3. Added 2  $\mu$ L pRSET to the tube
4. Flicked tube to mix and briefly centrifuged the samples down
5. Placed tube on ice for 30 minutes
6. Heat shocked E. coli in a 42 °C water bath for 30 seconds
7. Placed tubes on ice for 5 minutes
8. Added 475  $\mu$ L outgrowth medium to each tube
9. Incubated E. coli at 37 °C and shook at 400 rpm for 1 hour

**Results:** N/A

**Conclusion:** N/A

## **Plating G1 Transformed E. coli**

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Wednesday, May 31, 2023

**Protocol:**

1. Warmed an LB-agar ampicillin plate in a 37 °C incubator
2. Pipetted 50 µL E. coli cells onto each plate
3. Spread cells across plate using 3 glass beads for each one
4. Incubated at 37 °C overnight

**Results:** N/A

**Conclusion:** The colonies growing on the plate contain the pRSET vector because they were able to survive in the presence of ampicillin, and the pRSET vector confers ampicillin resistance.