

## Colony PCR

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Tuesday, June 27, 2023

**Protocol:**

**Find GoTaq protocol??**

1. Prepared a 163.2  $\mu\text{L}$  colony PCR reaction master mix with the following components:

Component	Volume
GoTaq <b>concentration?</b>	60 $\mu\text{L}$
Forward primer	1.6 $\mu\text{L}$
Reverse primer	1.6 $\mu\text{L}$
Nuclease-free water	80 $\mu\text{L}$
<b>Total</b>	<b>163.2<math>\mu\text{L}</math></b>

2. Pipetted 20  $\mu\text{L}$  master mix into 8 PCR tubes
3. Identified 8 colonies on the pRH plate circled each one with Sharpie
4. Swabbed each colony with a different sterile toothpick and swirled it in their respective PCR tube
5. Amplified the DNA in each 20  $\mu\text{L}$  reaction with a thermocycler using the following conditions:

Step	Temperature	Time	Number of cycles
Initial denaturation	98 °C	2 min	<b>1</b>
Denaturation	98 °C	30 sec	<b>30</b>
Annealing	57 °C	30 sec	<b>30</b>
Extension	72 °C	39 sec	<b>30</b>
Final extension	72 °C	2 min	<b>1</b>
Hold	10 °C	Forever	<b>-</b>

**Results:** N/A

**Conclusion:** N/A

## **Colony PCR Gel Purification (what exactly are you purifying?)**

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Tuesday, June 27, 2023

### **Protocol:**

1. Prepared a 1% agarose running gel:
  - a. Added 0.25 g of agarose to 25 mL 1X TBE in an Erlenmeyer flask
  - b. Heated flask in a microwave in 10 second increments until fully dissolved
  - c. Cooled flask until it was manageable to touch
  - d. Added 3  $\mu$ L SYBR Safe and swirled the flask to mix
  - e. Poured solution into a gel mold with a comb and left to solidify for 20 minutes
  - f. Removed comb and covered gel cast completely with 1X TBE
2. Loaded gel:
  - a. Loaded **4  $\mu$ L of 100 bp DNA ladder** into the first well
  - b. Loaded **20  $\mu$ L** colony PCR product into the second well
3. Placed gel cast in a blueGel electrophoresis system and left the sample to run for 30 minutes
4. **15 minutes after loading the 100 bp DNA ladder, a 1kb ladder was loaded into the same well**

### **Results:**

**Insert scanned photo (page 31)**

### **Conclusion:**

**3 and 4 have faint bands**

**7& 8 have bright bands – means that the colony PCR was a success – the bands mean that our primers landed inside of the bacteria which means that the insert is in the bacteria somewhere – primers should only bind to the insert.**

## **E. coli Inoculation**

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Tuesday, June 27, 2023

### **Protocol:**

1. Prepared 4 liquid cultures each with 3  $\mu$ L (are you sure its micro?) LB-ampicillin media
2. Swabbed colonies 3, 4, 7, and 8 with a different sterile toothpick and dropped it into their respective liquid culture
3. Incubated all 4 cultures at 37 °C and shook at 200 rpm overnight (18-20 hours?)

**Results:** N/A If there's a control, there should be no bacteria in the culture tube. This is quality assurance of the culture. The culture was not contaminated? Was there growth in the liquid cultures & checked May 3

**Conclusion:** N/A

Did you turn on a Bunsen burner while performing the inoculation?

Did you re flame the bottles and caps?

Supposed to re flame the culture tube and cap as well

Was there a negative control?