LB Broth Preparation

<u>Participants</u>: Tahmid Anwar, Pranav Bhavaraju, Rori Hoover, Patrick Jiang, Isabella (Bella) Lirtzman

Date: Tuesday, June 6, 2023

Protocol:

- 1. Autoclaved a 1 L glass bottle
- 2. Added 20 g Luria-Bertani broth powder to the bottle
- 3. Dissolved powder in 800 mL DI water
- 4. Loosely capped the glass bottle and autoclaved
- 5. Stored LB broth at 4°C

Results: N/A

Conclusion: N/A