

E. coli G1 DNA Miniprep

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Wednesday, May 3, 2023

Protocol:

Note: these steps follow those outlined in the Monarch® Plasmid DNA Miniprep Kit Protocol (NEB #T1010) General Guidelines by New England Biolabs.

1. Transferred 5 mL (what volume did you do, the notebook says 1-5 mL) of each liquid culture into separate microfuge tubes
2. Centrifuged all 4 tubes at 13,000 rpm for 30 seconds and discarded supernatant
3. Resuspended pelleted cells in each tube in 200 μ L Plasmid Resuspension Buffer and vortexed until fully resuspended
4. Added 200 μ L Plasmid Lysis Buffer to each tube and gently inverted the tubes until the solution turned dark pink
5. Incubated cells on the bench for 1 minute
6. Added 400 μ L Neutralization Buffer to each tube and inverted the tubes until the solution turned yellow and a precipitate formed
7. Incubated cells on the bench for 2 minutes
8. Centrifuged all 4 tubes at 13,000 rpm for 5 minutes (how many minutes did you centrifuge it? In the notebook you just wrote 2-5 minutes)
9. Inserted a spin column into 4 collection tubes and transferred the supernatant into their respective column
10. Centrifuged all 4 tubes at 13,000 rpm for 1 minute and discarded flow-through
11. Added 200 μ L Plasmid Wash Buffer 1 to each column and centrifuged all of the columns at 13,000 rpm for 1 minute
12. Added 400 μ L Plasmid Wash Buffer 2 to each column and centrifuged all of the columns at 13,000 rpm for 1 minute
13. Transferred each column to a clean 1.5 mL microfuge tube
14. Added 30 μ L Elution Buffer to each tube
15. Centrifuged all of the tubes at 13,000 rpm for 1 minute

Results: N/A

Conclusion: N/A

Measurement of **post-HiFi is G4** DNA Concentration

Participants: Brianna Branson, Rori Hoover, Patrick Jiang, Niam LeStourgeon

Date: Monday, June 5, 2023

Protocol:

1. Cleaned NanoDrop spectrophotometer with DI water
2. Blanked NanoDrop with 2 μL Elution Buffer
3. Cleaned NanoDrop before loading 2 μL sample
4. Repeated step 3 for all pRSET samples

Results:

Colony	DNA Concentration	260/280	260/230	Volume	Mass
3	243 ng/ μL	1.83	1.72	27.3 μL	6633.9 ng
4	252.6 ng/ μL	1.87	2.25	24.4 μL	6163.44 ng
7	234.3 ng/ μL	1.82	1.85	26.4 μL	6185.52 ng
8	236.0 ng/ μL	1.87	2.26	24.4 μL	5758.4 ng

Conclusion: N/A