

Re-Digest

- We're trying to see ~50 bp so we need ~1 µg DNA

Checklist

- G3.1: 13 µL
- G3.2: 15 µL
- G3.3: 15 µL
- G3.4: 13 µL
- G3.5: 13 µL
- G3.6: 15 µL
- G3.7: 18 µL
- G3.8: 24.5 µL

G3 Digestion

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Friday, July 28, 2023

Protocol:

1. Prepared 56 μL digestion master mix with the following components:

Component	Volume
10X rCutSmart buffer	40 μL
XbaI concentration?	8 μL
BamHI-HF	8 μL
Total	56 μL

2. Added 7 μL digestion master mix to each miniprep G3 sample
3. Prepared a 2% agarose running gel:
 - a. Added 0.4 g of agarose to 20 mL 1X TBE in an Erlenmeyer flask
 - b. Heated flask in a microwave in 10 second increments until fully dissolved
 - c. Cooled flask until it was manageable to touch
 - d. Added 3 μL SYBR Safe and swirled the flask to mix
 - e. Poured solution into a gel mold with a comb and left to solidify for 20 minutes
 - f. Removed comb and covered gel cast completely with 1X TBE
2. Loaded gel:
 - a. Loaded 3 μL 100 bp DNA ladder into the first well
 - b. Added 7 μL gel loading dye to each sample
 - c. Loaded ?? μL of each sample in ascending order from left to right starting with the second well
3. Placed gel cast in a blueGel electrophoresis system and left the sample to run for 1 hour

Results:

Insert scanned image from page 52

Conclusion: the PCR to remove His-tags did not work because