

E. coli G1 DNA Miniprep

Participants: Rori Hoover, Patrick Jiang

Date: Saturday, May 20, 2023

Protocol:

1. Centrifuged all 8 inoculated cultures at 13,000 rpm for 30 seconds and discarded supernatant
2. Resuspended pelleted cells in each tube in 200 μ L Plasmid Resuspension Buffer
3. Added 200 μ L Plasmid Lysis Buffer to each tube and gently inverted the tubes until the solution turned dark pink
4. Incubated cells on the bench for 1 minute
5. Added 400 μ L Neutralization Buffer to each tube and inverted the tubes until the solution turned yellow and a precipitate formed
6. Incubated cells on the bench for 2 minutes
7. Centrifuged all 8 tubes at 13,000 rpm for 5 minutes
8. Inserted a spin column into 8 collection tubes and transferred the supernatant into their respective column
9. Centrifuged all 8 tubes at 13,000 rpm for 1 minute and discarded flow-through
10. Added 200 μ L Plasmid Wash Buffer 1 to each column and centrifuged all of the columns at 13,000 rpm for 1 minute
11. Added 400 μ L Plasmid Wash Buffer 2 to each column and centrifuged all of the columns at 13,000 rpm for 2 minutes
12. Transferred each column to a clean 1.5 mL microfuge tube
13. Added 30 μ L Elution Buffer to each tube
14. Waited 1 minute before centrifuging all of the tubes at 13,000 rpm for 1 minute

Results: N/A

Conclusion: N/A