

Continuation of pRH DNA Purification

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Friday, June 23, 2023

Protocol:

Did you follow a company's protocol?

1. Centrifuged sample at # for 10 minutes at 0°C and removed supernatant
2. Added # μL 70% ethanol to the sample
3. Centrifuged at # for 2 minutes at 4°C and removed supernatant
4. Incubated sample on bench top with the tube open **for how long?**
5. Dissolved DNA pellet with 10 μL Elution Buffer

Results: N/A

Conclusion: N/A

Measurement of pRH DNA Concentration

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Monday, June 26, 2023

Protocol:

1. Cleaned NanoDrop spectrophotometer with DI water
2. Blanked NanoDrop with 2 μ L Elution Buffer
3. Cleaned NanoDrop before loading 2 μ L pRH DNA

Results: The pRH DNA concentration was found to be 17.6 ng/ μ L. The 260/280 ratio was 1.85, and the 260/230 ratio was 1.44.

Conclusion: **This is fine? Just a very low concentration**

G4.1, G4.2, and G4.3 HiFi DNA Assembly

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Monday, June 26, 2023

Protocol:

Note: these steps follow those outlined in the Barrick Lab protocol???

1. Prepared a 20 μL HiFi assembly reaction for each of the three inserts: G4.1, G4.2, and G4.3 with the following components:

| Component | Volume |
|---|------------------------------------|
| 0.056 pmols insert DNA (G4.1, G4.2, G4.3) | 5 μL |
| NEBuilder HiFi DNA Assembly Master Mix | 10 μL |
| Nuclease-free water | 5 μL |
| Total | 20 μL |

2. Prepared a 20 μL control reaction with the following components:

| Component | Volume |
|---|------------------------------------|
| 0.022 pmols pRSET (is the backbone pRSET or pRH?) | 5 μL |
| NEBuilder HiFi DNA Assembly Master Mix | 10 μL |
| Nuclease-free water | 5 μL |
| Total | 20 μL |

3. Did you amplify the DNA? What does 60 min 50°C mean?

| Step | Temperature | Time | Number of cycles |
|----------------------|-------------|---------|------------------|
| Initial denaturation | 98 °C | 30 sec | 1 |
| Denaturation | 98 °C | 10 sec | 30 |
| Annealing | 60 °C | 30 sec | 30 |
| Extension | 72 °C | 12 sec | 30 |
| Final extension | 72 °C | 2 min | 1 |
| Hold | 10 °C | Forever | - |

Results: N/A

Conclusion: N/A

pRH Gel Purification

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Monday, June 26, 2023

Protocol:

1. Prepared a 1% agarose running gel:
 - a. Added 0.2 g of agarose to 20 mL 1X TBE in an Erlenmeyer flask
 - b. Heated flask in a microwave in 10 second increments until fully dissolved
 - c. Cooled flask until it was manageable to touch
 - d. Added 2.5 μ L SYBR Safe and swirled the flask to mix
 - e. Poured solution into a gel mold with a comb and left to solidify for 20 minutes
 - f. Removed comb and covered gel cast completely with 1X TBE
2. Loaded gel:
 - a. Loaded 10 μ L of 1 kb DNA ladder into the first well
 - b. Added 2 μ L gel loading dye to the digested pRH
 - c. Loaded 2 μ L of pRH DNA into the second well
3. Placed gel cast in a blueGel electrophoresis system and left the sample to run for 30 minutes

Results:

Insert scanned photo (page 29)

Conclusion:

Inconclusive?

Took HiFi mix and ran it through a gel – ran too little of the HiFi and couldn't visualize anything

No follow up. Transformation does the same thing

Transformation of E. coli with **HiFi**

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Monday, June 26, 2023

Protocol:

1. Thawed NEB 10-beta Competent E. coli cells on ice for 10 minutes
2. Added 25 μ L E. coli cells to two microfuge tubes
3. Added 2 μ L pRH DNA??
4. Flicked tube to mix and briefly centrifuged the samples down
5. Placed tubes on ice for 30 minutes
6. Heat shocked E. coli in a 42 °C water bath for 30 seconds
7. Placed tubes on ice for 5 minutes
8. Added 475 μ L outgrowth medium to each tube
9. Incubated E. coli at 37 °C and shook at 400 rpm for 1 hour

Results: N/A

Conclusion: N/A

Plating ?? Transformed E. coli

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Monday, June 26, 2023

Protocol:

1. Warmed a LB-agar ampicillin plate in a 37 °C incubator
2. Pipetted 100 µL E. coli cells onto the plate
3. Spread cells across plate using 3 glass beads
4. Incubated at 37 °C overnight

Results: ???

Are there pictures??

Conclusion: ???

We did not transform a control so the danger is that gel purification is not perfect and contaminant whole plasmid has a much higher transfection efficiency that colonized things. We need to do colony PCR to determine if our insert is present.