

Gel

After discussion with Santiago, we realized that the error likely runs in the MegaWHOP and Transform. We are running a gel for digested and PCR products. Expect bands at ~7kb for the GZ, and ~6kb for the control PCR. The digested control expected no growth due to no cloning occurring in thermocycler, so all DNA was destroyed.

	1	2	3 (G1)	4
Order on gel:	NTL (Digested)	(GZ) (Digested)	NTL - GZ	
expected:	0 kb	~7kb	~6kb	~7kb
actual:	—	—	—	—

What is this?

mScarlet Megaprimer Generation

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Tuesday, May 16, 2023

Protocol:

1. Prepared a 50 μL reaction in a PCR tube with the following components:

Component	Volume
10 μM forward primer	2.5 μL
10 μM reverse primer	2.5 μL
10 ng template DNA (R8)	0.625 μL
Q5 High-Fidelity 2X Master Mix	25 μL
Nuclease-free water	19.4 μL
Total	50 μL

2. Amplified the DNA in the 25 μL reaction with a thermocycler, using the following conditions:

Step	Temperature	Time	Number of cycles
Initial denaturation	98 $^{\circ}\text{C}$	30 sec	1
Denaturation	98 $^{\circ}\text{C}$	10 sec	30
Annealing	61 $^{\circ}\text{C}$	30 sec	30
Extension	72 $^{\circ}\text{C}$	30 sec	30
Final extension	72 $^{\circ}\text{C}$	2 min	1
Hold	10 $^{\circ}\text{C}$	Forever	-

3. Prepared a 1% agarose running gel:


- a. Added 0.25 g of agarose to 25 mL 1X TBE in an Erlenmeyer flask
- b. Heated the flask in a microwave in 10 second increments until fully dissolved
- c. Cooled the flask until it was manageable to touch
- d. Added 2.5 μL SYBR Safe and swirled the flask to mix
- e. Poured the solution into a gel mold with a comb and left it to solidify for 20 minutes
- f. Removed the comb and covered the gel cast completely with 1X TBE

4. Loaded gel:

- a. Loaded 5 μL of 1 kb DNA ladder into the first well
- b. Added 3 μL gel loading dye to the PCR product
- c. Loaded 5 μL of PCR sample into the second well

5. Placed the gel cast in a blueGel electrophoresis system and left the sample to run for 30 minutes

Results:



Conclusion: We were successful in generating the mScarlet Megaprimer because a thin band appeared at 1 kb. At this point, we had still not realized that we had forgotten to dilute the forward and reverse primer from 100 μ M to 10 μ M, so the larger bands appearing at approximately 200 bp are again the result of self-primerization.

Gel Extraction of mScarlet Megaprimer

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Tuesday, May 16, 2023

Protocol:

1. Excised mScarlet Megaprimer DNA from the agarose gel and placed it in a microfuge tube.
2. Added 584.8 μL ethanol to the agarose gel piece
3. Incubated the agarose piece at 50 °C while shaking at 400 rpm until the gel completely dissolved
4. Inserted a spin column into a collection tube and loaded the sample into the column
5. Centrifuged the tube at 13,000 rpm for 1 minute and discarded the flow-through
6. Added 200 μL DNA Wash Buffer 1 to the column and centrifuged the column at 13,000 rpm for 1 minute
7. Added another 200 μL DNA Wash Buffer 1 to the column and centrifuged the column at 13,000 rpm for 1 minute
8. Transferred the spin column to a clean microfuge tube
9. Added 30 μL DNA Elution Buffer
10. Waited 1 minute before centrifuging the tube at 13,000 rpm for 1 minute

Results: N/A

Conclusion: N/A

Measurement of mScarlet Megaprimer DNA Concentration

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Tuesday, May 16, 2023

Protocol:

1. Cleaned NanoDrop spectrophotometer with DI water
2. Blanked NanoDrop with 2 μL Elution Buffer
3. Cleaned NanoDrop before loading 2 μL extracted mScarlet megaprimer

Results: The DNA concentration of the extracted mScarlet megaprimer sample was found to be 69.1 ng/ μL .

Conclusion: N/A

G1 with mScarlet Megaprimer Generation

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Tuesday, May 16, 2023

Protocol:

1. Prepared a 50 μ L MEGAWHOP reaction in a PCR tube with the following components:

Component	Volume
13.82 ng G1 plasmid	2 μ L
138.2 ng mScarlet megaprimer	1.56 μ L
Q5 High-Fidelity 2X Master Mix	25 μ L
Nuclease-free water	21.44 μ L
Total	50 μL

2. Amplified the DNA with a thermocycler using the following conditions:

Step	Temperature	Time	Number of cycles
Initial denaturation	98 °C	30 sec	1
Denaturation	98 °C	10 sec	30
Annealing	62 °C	30 sec	30
Extension	72 °C	3 min 35 sec	30
Final extension	72 °C	2 min	1
Hold	10 °C	Forever	-

Results: N/A

Conclusion: N/A