

E. coli Inoculation

Participants: Jennifer Brodsky (mentor)

Date: Tuesday, May 2, 2023

Protocol:

1. Prepared 8 liquid cultures each with ???
2. 8 colonies on the G1 plate were identified and circled in Sharpie
3. Added 3 mL kanamycin LB broth to a tube??
4. Swabbed each colony with a toothpick and added to their respective liquid culture
5. Incubated all 8 cultures at 37 °C and shook at 200 rpm overnight

Results: N/A

Conclusion: N/A