

E. coli pRSET-iGluSnFR3 DNA Miniprep

Participants: Rori Hoover, Patrick Jiang

Date: Saturday, June 3, 2023

Protocol:

Note: these steps follow those outlined in the Monarch® Plasmid DNA Miniprep Kit Protocol (NEB #T1010) General Guidelines by New England Biolabs.

1. Transferred 5 mL of each liquid culture into 8 different microfuge tubes
2. Centrifuged all 8 tubes at 13,000 rpm for 30 seconds and discarded supernatant
3. Resuspended pelleted cells in each tube in 200 μ L Plasmid Resuspension Buffer
4. Added 200 μ L Plasmid Lysis Buffer to each tube and gently inverted the tubes until the solution turned dark pink
5. Incubated cells on the bench for 1 minute
6. Added 400 μ L Neutralization Buffer to each tube and inverted the tubes until the solution turned yellow and a precipitate formed
7. Incubated cells on the bench for 2 minutes
8. Centrifuged all 8 tubes at 13,000 rpm for 5 minutes
9. Inserted a spin column into 8 collection tubes and transferred the supernatant into their respective column
10. Centrifuged all 8 tubes at 13,000 rpm for 1 minute and discarded flow-through
11. Added 200 μ L Plasmid Wash Buffer 1 to each column and centrifuged all of the columns at 13,000 rpm for 1 minute
12. Added 400 μ L Plasmid Wash Buffer 2 to each column and centrifuged all of the columns at 13,000 rpm for 2 minutes
13. Transferred each column to a clean 1.5 mL microfuge tube
14. Added 30 μ L Elution Buffer to each tube
15. Waited 1 minute before centrifuging all of the tubes at 13,000 rpm for 1 minute

Results: N/A

Conclusion: We will measure the DNA concentration of each miniprep sample for downstream use.