

## **Transformation of E. coli with G2-D**

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Friday, June 30, 2023

**Protocol:**

1. Thawed NEB 10-beta Competent E. coli cells on ice for 10 minutes
2. Added 25  $\mu\text{L}$  E. coli cells to a microfuge tube
3. Added 2  $\mu\text{L}$  G2-D plasmid to the tube
4. Flicked tube to mix and briefly centrifuged the samples down
5. Placed tubes on ice for 30 minutes
6. Heat shocked E. coli in a 42 °C water bath for 30 seconds
7. Placed tubes on ice for 5 minutes
8. Added 475  $\mu\text{L}$  outgrowth medium to each tube
9. Incubated E. coli at 37 °C and shook at 400 rpm for 1 hour

**Results:** N/A

**Conclusion:** N/A

**No control?**

**Chemistry team – as practice**

## **Plating G2-D Transformed E. coli**

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Friday, June 30, 2023

**Protocol:**

1. Warmed 2 LB-agar kanamycin plates in a 37 °C incubator
2. Pipetted 100 µL E. coli cells onto a plate
3. Spread cells across plate using 3 glass beads
4. Incubated at 37 °C overnight

**Results:** There was an overabundance of pink growth on the plate with E. coli that had G2-D transformed into them.

**Conclusion:** Our results indicate that we were successful in transforming E. coli.