

## G3 Plasmid Digestion

**Participants:** Dr. Benjamin Bartelle (mentor)

**Date:** Saturday, July 15, 2023

**Protocol:**

1. Added 0.5  $\mu$ L DpnI to each PCR tube
2. Incubated samples for 1 hour at room temperature

**Results:** N/A

**Conclusion:** N/A

DpnI destroys methylated DNA – should not eat synthesized DNA from PCR

## Transformation of E. coli with G3

**Participants:** Dr. Benjamin Bartelle (mentor)

**Date:** Saturday, July 15, 2023

**Protocol:**

1. Thawed what kind? NEB E. coli cells on ice for 10 minutes
2. Added 25  $\mu\text{L}$  E. coli cells to a microfuge tube
3. Added 2  $\mu\text{L}$  G3 to the tube
4. Flicked tube to mix and briefly centrifuged the samples down
5. Placed tubes on ice for 30 minutes
6. Heat shocked E. coli in a 42 °C water bath for 30 seconds
7. Placed tubes on ice for 5 minutes
8. Added 475  $\mu\text{L}$  outgrowth medium to each tube
9. Incubated E. coli at 37 °C and shook at 400 rpm for 1 hour

**Results:** N/A

**Conclusion:** N/A

Was there a control? Idk if this protocol is correct. If it's BL21, it was done differently by Bella and Brianna

## **Plating G3 Transformed E. coli**

**Participants:** Dr. Benjamin Bartelle (mentor)

**Date:** Saturday, July 15, 2023

**Protocol:**

1. Warmed an LB-agar kanamycin plate in a 37 °C incubator
2. Serially streaked G3 transformed E. coli onto the plate using a sterilized streaking tool
3. Incubated at 37 °C over the weekend

**Results:** ??

**Conclusion:** Growth on the plate