

## Colony PCR

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Monday, July 24, 2023

**Protocol:**

**Did you follow a company's protocol?**

1. Prepared 2 29.5  $\mu\text{L}$  colony PCR reactions with the following components:

Component	Volume
GoTaq <b>concentration?</b>	20 $\mu\text{L}$
Forward primer	0.25 $\mu\text{L}$
Reverse primer	0.25 $\mu\text{L}$
Nuclease-free water	9 $\mu\text{L}$
<b>Total</b>	<b>29.5 <math>\mu\text{L}</math></b>

2. Amplified the DNA in each 29.5  $\mu\text{L}$  reaction with a thermocycler using the following conditions:

Step	Temperature	Time	Number of cycles
Initial denaturation	98 °C	2 min	1
Denaturation	98 °C	30 sec	30
Annealing	57 °C	30 sec	30
Extension	72 °C	39 sec	30
Final extension	72 °C	2 min	1
Hold	10 °C	Forever	-

**Results:** N/A

**Conclusion:** N/A

## **Colony PCR Gel Purification (what exactly are you purifying?)**

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Monday, July 24, 2023

### **Protocol:**

1. Prepared a 1% agarose running gel:
  - a. Added 0.2 g of agarose to 20 mL 1X TBE in an Erlenmeyer flask
  - b. Heated flask in a microwave in 10 second increments until fully dissolved
  - c. Cooled flask until it was manageable to touch
  - d. Added 3  $\mu$ L SYBR Safe and swirled the flask to mix
  - e. Poured solution into a gel mold with a comb and left to solidify for 20 minutes
  - f. Removed comb and covered gel cast completely with 1X TBE
2. Loaded gel:
  - a. Loaded #  $\mu$ L of 100 bp DNA ladder into the first well
  - b. Added #  $\mu$ L gel loading dye to colony PCR product
  - c. Loaded #  $\mu$ L colony PCR product into the second well
3. Placed gel cast in a blueGel electrophoresis system and left the sample to run for 30 minutes
4. 15 minutes after loading the 100 bp DNA ladder, a 1kb ladder was loaded into the same well

### **Results:**

**Insert scanned photo (page 40)**

### **Conclusion:**

**Colony PCR success on small colony – insert was unclear? – what is your conclusion?**

### **G3 Assembly**

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Monday, July 24, 2023

**Protocol:**

**5  $\mu$ L of PCR (add DpnI)**

**Taq**

**PCR: 5  $\mu$ L**

**Master mix: 5  $\mu$ L**

**10  $\mu$ L total**

**This is a test if the annealing temperature was an issue**

**Results:** ??

**Conclusion:** ??

## **Transformation of E. coli with G3**

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Monday, July 24, 2023

**Protocol:**

1. Thawed **what kind?** E. coli cells on ice for 10 minutes
2. Added 25  $\mu\text{L}$  E. coli cells to a microfuge tube
3. Added 2  $\mu\text{L}$  G3 to the tube
4. Flicked tube to mix and briefly centrifuged the samples down
5. Placed tubes on ice for 30 minutes
6. Heat shocked E. coli in a 42 °C water bath for 30 seconds
7. Placed tubes on ice for 5 minutes
8. Added 475  $\mu\text{L}$  outgrowth medium to each tube
9. Incubated E. coli at 37 °C and shook at 250 rpm for 1 hour

**Results:** N/A

**Conclusion:** N/A

Was there a control? Idk if this protocol is correct. If it's BL21, it was done differently by Bella and Brianna

## Plating G3 Transformed E. coli

**Participants:** Brianna Branson, Rori Hoover, Patrick Jiang

**Date:** Monday, July 24, 2023

**Protocol:**

1. Warmed an **LB-agar kanamycin** plate in a 37 °C incubator
2. Serially streaked G3 transformed E. coli onto the plate using a sterilized streaking tool
3. Incubated at 37 °C over the weekend
- 4. Plate 50 µL??**

**Results:** **Nothing grew on the G3 plate**

**Conclusion:** **To troubleshoot, we will spin down the unplated transformed cells, discard the supernatant, resuspend with 100 µL outgrowth media, and plate 100 µL**

**Why did you decide to do this to troubleshoot**