

## **Sample Preparation for Sequencing**

**Participants:** Rori Hoover, Patrick Jiang

**Date:** Monday, May 22, 2023

**Protocol:**

1. Added 3.28  $\mu\text{L}$  91.4 ng/ $\mu\text{L}$  plasmid (Colony D) to a microfuge tube
2. Diluted plasmid with 6.82  $\mu\text{L}$  nuclease-free water
3. Shipped sample to Plasmidosaurus for sequencing

**Results:** N/A

**Conclusion:** N/A