G5 Amplification??

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Friday, July 28, 2023

Protocol:

1. Prepared a 25 μ L reaction in a PCR tube with the following components:

Component	Volume	
10 μM forward primer	1.25 μL	
10 μM reverse primer	1.25 μL	
10 ng template DNA (??)	1 μL	
Q5 High-Fidelity 2X Master Mix	12.5 μL	
Nuclease-free water	<mark>4 μL</mark>	
Total	25 μL	

2. Amplified the DNA in the 25 μL reaction with a thermocycler using the following conditions:

Step	Temperature	Time	Number of cycles
Initial denaturation	98 °C	30 sec	1
Denaturation	98 °C	10 sec	30
Annealing	<mark>57</mark> °C	30 sec	30
Extension	72 °C	30 sec	30
Final extension	72 °C	2 min	1
Hold	10 °C	Forever	-

Results: N/A

Conclusion: N/A

Plating on M9 with Hexanoic Acid in Two Forms

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Friday, July 28, 2023

Protocol:

- 1. Warmed M9 plates in a 37 °C incubator
- 2. Serially streaked G4C?? transformed E. coli onto the plate using a sterilized streaking tool
- 3. Incubated at 37 °C over the weekend

Results: ???

Conclusion: ???