G4.4 E. coli Inoculation

Participants: Brianna Branson, Rori Hoover, Patrick Jiang

Date: Tuesday, July 11, 2023

Protocol:

- 1. Prepared 3 liquid cultures with ampicillin media?? using a serological pipette
- 2. Identified 3 colonies on the G4.4 plate and circled each one with Sharpie
- 3. Swabbed each colony with a different sterile toothpick and dropped it into their respective liquid culture
- 4. Incubated all 3 cultures at 37 °C and shook at 200 rpm overnight

Results: N/A If there's a control, there should be no bacteria in the culture tube. This is quality assurance of the culture. The culture was not contaminated? Was there growth in the liquid cultures & checked May 3

Conclusion: N/A

Did you turn on a Bunsen burner while performing the inoculation?

Did you reflame the bottles and caps?

Supposed to reflame the culture tube and cap as well

Was there a negative control?